



I-STAT Training - Quick Tips

Point-Of-Care-Testing Department
Clinical Laboratories
Stanford Health Care
Stanford Children's Health



Patient/ Operator ID entry

Patient ID entry:

- Joint Commission requirement for 2 forms Patient Identification.
- Scan Patient armband for CSN = 9-digit in LPCH and 12-digit in SHC
- Enter twice if manually enter in.
- **Operator/User ID:**
 - Trained personnel will be issued a bar-coded ID label for scanning with the ISTAT meter.
 - It is the 5-digit or 6-digit of the ID badge number.



I-STAT cartridge storage, handling and supplies

- **Cartridges are stored in refrigeration at 2 - 8° C (35 - 46° F). Use within manufacturer expiration date listed on the box. It takes 5 minutes to warm up one packet and one hour to warm up one whole box of cartridges.**
- **Different cartridges have different expiration date at room temp storage. Write room temp expiration date on individual pouch or box if pouches are kept inside the box.**
- **Handle cartridges by the edges only. Burst or damaged cartridge will not give results.**

Your unit clerk has to go to the clinical lab for cartridge supplies - call ahead to **650-736-2127
Blood Gas Lab.**



Blood sample information

- **Non-Heparinized 1CC TB syringe. Blood has to be tested within 3 minutes.**
- **Heparinized syringe. Use within 10 minutes with blood gently mixed before sampling.**
- **For coagulation tests, i.e. ACT, test the sample immediately.**
- **Do not use blood if it is clotted**
- **Do not draw air bubbles into the syringe – blood gas values are affected.**
- **Perform test promptly**

All of the above factors affect blood gas results.



I-STAT testing process

- It takes 2 minutes to see results. Coagulation testing take up to 10 minutes to complete.
- Do not remove cartridge during testing process (Cartridge Locked) – instrument will be damaged if cartridge is forcibly removed
- Instrument should maintain at a horizontal level during testing.
- Test process and error messages are shown on the bottom of the screen.
- Errors with *******, **<** or **>** will allow any results to transmit to Epic.
- Clean meter with Sani-Cloth wipe after each patient use per Infection Control protocol

Interferences for I-STAT testing

- Syringes with >10 units of Heparin falsely lowers ionized calcium results.
- Abnormally low or high Total Proteins in patient can falsely decrease or increase the Hematocrit value.
Hct values are used for trending only, not as an indicator for transfusion.
- Do not remove cartridge during testing process (**Cartridge Locked**) – instrument will be damaged if cartridge is forcibly removed.


When testing needs to be repeated

- **Conditions that a testing needs to be repeated with a fresh blood sample:**
 - **Sample is clotted**
 - **Sample is contaminated with IV fluids**
 - **Abnormal results that do not match with patient treatment protocol. Alert values must be documented and reported to the attending physicians to ensure proper patient treatment.**
 - **Results with ***, >, < or <> - if persistent do not use ISTAT for testing. Send sample to the Lab. These errors are transmitted to EPIC system.**

If tests are repeated within 5 minutes, only the second set of results are allowed to stay as the permanent record.

Erroneous results can be deleted by pressing “ 8 ” in Field 2 under Sample Type screen or call POCT at 4-8934.

Do not use i-STAT for patient testing if repeated *, < or > are shown. Send blood samples to the Clinical Lab for analysis.**



Quality Controls (QC) for I-STAT: Internal, external and liquid controls

- **Internal Simulator (QC):**
 - It is programmed to run automatically every 8 hours of instrument use.
- **External Simulator (QC):**
 - Use the metallic QC when prompted by the instrument. Do not touch the metallic sensors.
 - Use this External QC to test meter is it is dropped to ensure instrument integrity. Patient testing is disable if meter is broken. Call POCT # 724-8934 to retrieve broken meter.
- **Liquid Controls:**
 - Liquid QC is performed in the units monthly
 - Liquid QC is also required on each new shipment/Lot.




Training and competency

- Each I-STAT user must have **Initial training and competency.**
- A **6-month competency** and then an **Annual competency** are required.
- **Hands-on demonstration** is required as part of the competencies.

Training and competency

- **When the competency certification expires within 30-day timeframe, a message appears on the display screen to alert you the expiration date.**



- **Users contact POCT 4-8934 Option 4 to complete the required test before the due date.**



Training and competency

- **User will be locked out for testing if competencies are not kept up to date.**
- **Re-training is required once the user ID expired in the ISTAT database.**
- **Sharing operator ID is not permissible as it may become a patient safety issue.**

Troubleshooting the I-STAT instrument

- Avoid test errors to reduce cost and to improve efficiency with patient care.
- I-STAT instrument performs quality checks for each testing.
- Examples of error messages and corrective actions:
 - “**Low Batteries or Dead Batteries**” – Replace the 9 volt batteries
 - “**External Simulator Run Required**” – Insert the metallic EQC to run quality checks
 - “**Analyzer Error**” – Place the instrument in the docking station to reset software profile. Run the electronic Simulator afterwards.
 - “**Cartridge Error**” – Try another cartridge, previous cartridge may be mishandled
 - “**Invalid or Expired CLEW**” – Call POCT department @ 650-724-8934 to update I-STAT software



Troubleshooting the I-STAT - Continue

- **“Unable to Position Sample”** – Snap closure has not been closed tightly or sample has clotted
 - **“Cartridge Pre-burst”** – Calibration pouch of the I-STAT cartridge has been damaged, use another cartridge
 - **“Sample underfilled or overfilled”** – Blood should be sampled up to the blue arrow point on the cartridge, use another cartridge
-
- **Notify POCT department 650-724-8934 Option 4 for I-STAT problems and replacements.**
 - **Notify SHC IT 3-3333 or LPCH IT 8-7500 if the results are not transmitted to EPIC.**
 - **Do not take ISTAT to Material Management or Clinical Engineering.**

Monthly Compliance Report – Joint Commission Requirements

Monthly Compliance checks:

- ✓ **User training/competency/certification status.**
- ✓ **Invalid patient ID when performing the i-STAT test.**
- ✓ **Meter(s) are cleaned per Infection Control protocol.**
- ✓ **Room and refrigerator temperatures are monitored and recorded daily on chart/log.**
- ✓ **ISTAT cartridges are properly **dated** when store at room temperature.**
- ✓ **ISTAT cartridges are properly **disposed** when expired with room temp storage or manufacturer expiration date.**
- ✓ **The acceptable monthly compliance score is >90%. Reports are sent to unit medical director.**